

Quick Reference Guide

ViewRNA™ ISH Tissue 1-Plex Assay

! **IMPORTANT:** If running ViewRNA ISH Tissue 1-Plex Assay for the first time, please refer to the *ViewRNA ISH Tissue 1-Plex Assay User Manual* to review assay guidelines and detailed procedures.

Part 1: Sample Preparation and Target Probe Hybridization

Step	Task
1 Bake Slides 65 min	<ul style="list-style-type: none"> ■ Label slides with pencil ■ Bake slides at 60 ± 1 °C for 60 min
2 Preparation for Part 1	<ul style="list-style-type: none"> ■ Verify hybridization system is set to 40 ± 1 °C and appropriately humidified ■ Prepare: <ul style="list-style-type: none"> □ 2 L 1X PBS □ 200 mL 10% NBF in 1X PBS □ 3 L Wash Buffer □ 500 mL 1X Pretreatment Solution □ 200 mL Storage Buffer (optional) □ 1 L ■ Ensure availability of: <ul style="list-style-type: none"> □ 400 mL 100% ethanol □ 1.4 L ddH₂O □ 600 mL xylene or 400 mL HistoClear □ 200 mL Gill's Hematoxylin I □ 200 mL of 3 µg/mL DAPI in 1X PBS (optional, for fluorescence detection) ■ Thaw probe set(s) and place on ice ■ Prewarm 40 mL 1X PBS and Probe Set Diluent QF to 40 °C <ul style="list-style-type: none"> ■ Optional for 1-day assay: <ul style="list-style-type: none"> □ Prewarm Amplifier Diluent QF and Label Probe Diluent QF to 40 °C □ Thaw PreAmp 1 QF and Amp1 QF, Mix, briefly spin down to collect contents, and place on ice until use □ Briefly spin down the Label Probe-AP and place on ice □ Bring Fast Red Tablets, Napthol Buffer, and AP Enhancer Solution to RT □ Prepare 1 L 0.01% ammonium hydroxide under a fume hood
3 Deparaffinization 30 min	<p>If using xylene (work in a fume hood):</p> <ul style="list-style-type: none"> ■ Incubate slides 3 times, each for 5 min in xylene ■ Wash slides 2 times in 100% ethanol, 5 min each wash ■ Decant 100% ethanol from slides and let air dry <p>If using Histo-Clear:</p> <ul style="list-style-type: none"> ■ Bake slides at 80 ± 1 °C for 3 min ■ Incubate slides 2 times, 5 min each in HistoClear ■ Wash slides 2 times in 100% ethanol, 5 min each wash ■ Decant 100% ethanol from slides and let air dry
4 Draw Hydrophobic Barrier 40 min	<ul style="list-style-type: none"> ■ Create hydrophobic barrier ■ Allow slides to air dry at RT for 20-30 min ■ Heat 1X Pretreatment Solution to 90-95 °C while slides are drying
5 Heat Pretreatment 10-25 min, depending on optimized time	<ul style="list-style-type: none"> ■ Heat slides in 1X Pretreatment Solution at 90-95 °C for the optimal time determined in the pretreatment assay optimization procedure ■ Wash slides 2 times in ddH₂O, 1 min each wash ■ Transfer slides to 1X PBS <p>IMPORTANT: Do not let the slides dry out from this point forward.</p>
6 Protease Digestion and Fixation 30-50 min, depending on optimized time	<ul style="list-style-type: none"> ■ Prepare 1:100 working protease solution in prewarmed 1X PBS ■ Add working protease solution to slides ■ Incubate at 40 °C for optimal time determined in the pretreatment assay optimization procedure ■ Wash slides 2 times in 1X PBS, 1 min each wash ■ Fix slides in 10% NBF at RT for 5 min under a fume hood ■ Wash slides 2 times in 1X PBS, 1 min each wash

Step	Task
7 Target Probe Set Hybridization 3 hr 10 min	<ul style="list-style-type: none"> ■ Prepare 1:50 working probe set solution in prewarmed Probe Set Diluent QF ■ Add working probe set solution to slides ■ Incubate at 40 °C for 3 hr
8 Wash Slides 8 min	Wash slides 3 times in Wash Buffer, 2 min each wash
9 Optional Stop Point 1 min	<ul style="list-style-type: none"> ■ Store slides in Storage Buffer at RT for up to 24 hr. Cover dish with lid or sealing film to prevent evaporation. ■ Store remaining 1X PBS and Wash Buffer at RT for use in Part 2

Part 2: Signal Amplification and Detection

Step	Task
10 Preparation for Part 2 10 min	<ul style="list-style-type: none"> ■ Pour Gill's Hematoxylin I into a clear staining dish, store at RT protected from light. ■ Optional – Prepare 200 mL 3 µg/mL DAPI, store at 4 °C until use or place on ice ■ Prewarm Amplifier Diluent QF and Label Probe Diluent QF to 40 °C ■ Thaw PreAmp1 QF and Amp1 QF, place on ice ■ Bring Fast Red Tablets, Naphthol Buffer, and AP Enhancer to RT ■ Prepare 1 L 0.01% ammonium hydroxide
11 Wash Slides 8 min	Remove slides from Storage Buffer and wash 3 times in Wash Buffer, 2 min each wash with constant, vigorous agitation
12 PreAmplifier Hybridization 35 min	<ul style="list-style-type: none"> ■ Prepare 1:100 Working PreAmp1 Solution in prewarmed Amplifier Diluent QF ■ Add working PreAmp1 solution to slides ■ Incubate at 40 °C for 25 min
13 Wash Slides 8 min	Wash slides 3 times in Wash Buffer, 2 min each wash with constant, vigorous agitation
14 Amplifier Hybridization 20 min	<ul style="list-style-type: none"> ■ Prepare 1:100 working Amp1 solution in prewarmed Amplifier Diluent QF ■ Add working Amp1 solution to slides ■ Incubate at 40 °C for 15 min
15 Wash Slides 8 min	Wash slides 3 times in Wash Buffer, 2 min each wash
16 Label Probe-AP Hybridization 20 min	<ul style="list-style-type: none"> ■ Prepare 1:1000 working Label Probe-AP solution in prewarmed Label Probe Diluent QF ■ Add working Label Probe-AP solution to slides ■ Incubate at 40 °C for 15 min
17 Wash Slides 12 min	Wash slides 3 times in Wash Buffer, 3 min each wash
18 Apply Fast Red Substrate 45 min	<ul style="list-style-type: none"> ■ Add AP-Enhancer to slides and incubate at RT for 5 min ■ Prepare Fast Red Substrate (1 Fast Red Tablet/5 mL Naphthol Buffer) ■ Decant AP-Enhancer and add Fast Red Substrate to slides, incubate at 40 °C for 30 min ■ Rinse slides in 1X PBS for 1 min
19 Counterstain 25 min	<ul style="list-style-type: none"> ■ Incubate slides in Gill's Hematoxylin I stain at RT for 5-10 sec ■ Wash slides 3 times in ddH₂O, 1 min each wash ■ Incubate in 0.01% ammonium hydroxide at RT for 10 sec ■ Wash slides in ddH₂O for 1 min ■ Optional – Incubate slides in DAPI at RT for 1 min, wash slides in ddH₂O for 1 min ■ Let slides completely air dry at RT (~20 min)

Step	Task
20 Mount and Image 40 min	<p>DAKO Ultramount:</p> <p>For no coverslipping</p> <ul style="list-style-type: none"> ■ Add Ultramount to tissue sections ■ Place slides in a 70 °C oven/incubator for 10-30 min ■ Observe under brightfield or fluorescence microscope ■ Store slides at RT <p>For post mounting with coverslip</p> <ul style="list-style-type: none"> ■ Work under a fume hood and follow the no coverslip procedure ■ Allow the slides to come to RT ■ Apply Histomount directly on top of the dried Ultramount ■ Place coverslip ■ Air dry at RT for 15 min ■ Observe under brightfield or fluorescence microscope ■ Store slides at RT <p>Innovex Advantage Mounting Media:</p> <ul style="list-style-type: none"> ■ Add Advantage Mounting Media to cover glass and invert tissue slide to cover ■ Flip over, allow slides to dry ■ Seal all four edges with nail polish ■ Observe under brightfield or fluorescence microscope ■ Store slides at RT

www.affymetrix.com/panomics

For research use only. Not for use in diagnostic procedures.