

Quick Reference Guide

ViewRNA™ ISH Tissue 2-Plex Assay

FFPE Samples

! **IMPORTANT:** If running ViewRNA ISH Tissue 2-Plex Assay for the first time, please refer to the *ViewRNA ISH Tissue 2-Plex Assay User Manual* to review assay guidelines and detailed procedures.

Part 1: Sample Preparation and Target Probe Hybridization

Step	Task
1 Bake Slides 65 min	<ul style="list-style-type: none"> ■ Label slides with pencil ■ Bake slides at 60 °C for 60 min
2 Preparation for Part 1	<ul style="list-style-type: none"> ■ Verify hybridization system is set to 40 °C and appropriately humidified ■ Prepare: <ul style="list-style-type: none"> □ 3 L 1X PBS □ 200 mL 10% NBF in 1X PBS □ 4 L Wash Buffer □ 500 mL 1X Pretreatment Solution □ 200 mL Storage Buffer ■ Ensure availability of: <ul style="list-style-type: none"> □ 600 mL 100% ethanol □ 1.4 L ddH₂O □ 600 mL xylene or 400 mL Histo-Clear □ 200 mL Gill's Hematoxylin I □ 200 mL of 3 µg/mL DAPI in 1X PBS (optional, for fluorescence detection) ■ Thaw probe set(s) on ice ■ Prewarm 40 mL 1X PBS and Probe Set Diluent QF to 40 °C <ul style="list-style-type: none"> ■ Optional for 1-day assay: <ul style="list-style-type: none"> □ Prewarm PreAmplifier Mix QT, Amplifier Mix QT, and Label Probe Diluent QF to 40 °C □ Briefly spin down Label Probe 1-AP, Label probe 6-AP, and Blue Reagents, place on ice □ Bring Fast Red Tablets, Naphthol Buffer, Blue Buffer, and AP Enhancer Solution to RT □ Prepare 1 L 0.01% ammonium hydroxide under a fume hood
3 Deparaffinization 30 min	<p>If using xylene (work in a fume hood):</p> <ul style="list-style-type: none"> ■ Incubate slides 3 times each for 5 min in xylene ■ Wash slides 2 times in 100% ethanol, 5 min each wash ■ Decant 100% ethanol from slides and let air dry <p>If using Histo-Clear:</p> <ul style="list-style-type: none"> ■ Bake slides at 80 ± 1 °C for 3 min ■ Incubate slides 2 times, 5 min each in HistoClear ■ Wash slides 2 times in 100% ethanol, 5 min each wash ■ Decant 100% ethanol from slides and let air dry
4 Draw Hydrophobic Barrier 40 min	<ul style="list-style-type: none"> ■ Create hydrophobic barrier ■ Allow slides to air dry at RT for 20-30 min ■ Heat 1X Pretreatment Solution to 90-95 °C while slides are drying
5 Heat Pretreatment 10-25 min, depending on optimized time	<ul style="list-style-type: none"> ■ Heat slides in 1X Pretreatment Solution at 90-95 °C for the optimal time determined in the pretreatment assay optimization procedure ■ Wash slides 2 times in ddH₂O, 1 min each wash ■ Transfer slides to 1X PBS <p>IMPORTANT: Do not let the slides dry out from this point forward.</p>
6 Protease Digestion and Fixation 30-50 min, depending on optimized time	<ul style="list-style-type: none"> ■ Prepare 1:100 working protease solution in prewarmed 1X PBS ■ Add working protease solution to slides ■ Incubate at 40 °C for optimal time determined in the pretreatment assay optimization procedure ■ Wash slides 2 times in 1X PBS, 1 min each wash ■ Fix slides in 10% NBF at RT for 5 min under a fume hood ■ Wash slides 2 times in 1X PBS, 1 min each wash

Step	Task
7 Target Probe Set Hybridization 2 hr 10 min	<ul style="list-style-type: none"> ■ Prepare 1:40 working probe set solution in prewarmed Probe Set Diluent QT ■ Add working probe set solution to slides ■ Incubate at 40 °C for 2 hr
8 Wash Slides 8 min	Wash slides 3 times in Wash Buffer, 2 min each wash
9 Optional Stop Point 1 min	<ul style="list-style-type: none"> ■ Store slides in Storage Buffer at RT for up to 24 hr. Cover dish with lid or sealing film to prevent evaporation. ■ Store 1X PBS and Wash Buffer at RT for use in Part 2

Part 2: Signal Amplification and Detection

Step	Task
10 Preparation for Part 2 10 min	<ul style="list-style-type: none"> ■ Pour Gill's Hematoxylin I into a clear staining dish, store at RT protected from light. ■ Optional – Prepare 200 mL 3 µg/mL DAPI, store at 4 °C until use or place on ice ■ Prewarm PreAmplifier Mix QT, Amplifier Mix QT and Label Probe Diluent QF to 40 °C ■ Place Label Probe 1-AP, Label Probe 6-AP, and Blue reagents on ice ■ Bring Fast Red Tablets, Naphthol Buffer, AP Enhancer, and Blue Buffer to RT ■ Prepare 1 L 0.01% ammonium hydroxide in ddH₂O
11 Wash Slides 8 min	Remove slides from Storage Buffer and wash 3 times with Wash Buffer with constant and vigorous agitation, 2 min each wash
12 PreAmplifier Hybridization 35 min	Add PreAmplifier Mix QT directly to slides, incubate at 40 °C for 25 min
13 Wash Slides 8 min	Wash slides 3 times in Wash Buffer with constant and vigorous agitation, 2 min each wash
14 Amplifier Hybridization 20 min	Add Amplifier Mix QT directly to slides, incubate at 40 °C for 15 min
15 Wash Slides 8 min	Wash slides 3 times in Wash Buffer with constant and vigorous agitation, 2 min each wash
16 Label Probe 6-AP Hybridization 20 min	<ul style="list-style-type: none"> ■ Prepare 1:1000 working Label Probe 6-AP solution ■ Add working Label Probe 6-AP solution to slides, incubate at 40 °C for 15 min
17 Wash Slides 12 min	Wash slides 3 times in Wash Buffer with constant and vigorous agitation, 3 min each wash
18 Apply Fast Blue Substrate 35 min	<ul style="list-style-type: none"> ■ Prepare Fast Blue Substrate ■ Add Fast Blue Substrate to slides, incubate at RT in the dark for 30 min
19 Wash Slides 12 min	Wash slides 3 times in Wash Buffer with constant and vigorous agitation, 3 min each wash
20 Quench Label Probe 6-AP 35 min	<ul style="list-style-type: none"> ■ Add AP Stop QT to slides, incubate at RT in the dark for 30 min ■ Wash slides 2 times in 1X PBS, 1 min each wash ■ Rinse slides in Wash Buffer for 1 min
21 Label Probe 1-AP Hybridization 20 min	<ul style="list-style-type: none"> ■ Prepare 1:1000 working Label Probe 1-AP solution in prewarmed Label Probe Diluent QF ■ Add working Label Probe 1-AP solution to slides, incubate at 40 °C for 15 min
22 Wash Slides 12 min	Wash slides 3 times in Wash Buffer, 3 min each wash

Step	Task
23 Apply Fast Red Substrate 45 min	<ul style="list-style-type: none"> ■ Add AP-Enhancer to slides and incubate at RT for 5 min ■ Prepare Fast Red Substrate (1 Fast Red Tablet/5 mL Naphthol Buffer) ■ Decant AP-Enhancer and add Fast Red Substrate to slides, incubate at 40 °C for 30 min ■ Wash slides in 1X PBS for 1 min
24 Counterstain 25 min	<ul style="list-style-type: none"> ■ Incubate slides in Gill's Hematoxylin I stain at RT for 5-10 sec ■ Wash slides 3 times in ddH₂O, 1 min each wash ■ Incubate in 0.01% ammonium hydroxide at RT for 10 sec ■ Wash slides in ddH₂O for 1 min ■ Optional – Incubate slides in DAPI at RT for 1 min, wash slides in ddH₂O for 1 min ■ Let slides completely air dry at RT (~20 min)
25 Mount and Image 40 min	<p>DAKO Ultramount: For no coverslipping</p> <ul style="list-style-type: none"> ■ Add Ultramount to tissue sections ■ Place slides in a 70 °C oven/incubator for 10-30 min ■ Observe under brightfield or fluorescence microscope ■ Store slides at RT <p>For post mounting with coverslip</p> <ul style="list-style-type: none"> ■ Work under a fume hood and follow the no coverslip procedure ■ Allow the slides to come to RT ■ Apply Histomount directly on top of the dried Ultramount ■ Place coverslip ■ Air dry at RT for 15 min ■ Observe under brightfield or fluorescence microscope ■ Store slides at RT <p>Innovex Advantage Mounting Media:</p> <ul style="list-style-type: none"> ■ Add Advantage Mounting Media to cover glass and invert tissue slide to cover ■ Flip over, allow slides to dry ■ Seal all four edges with nail polish ■ Observe under brightfield or fluorescence microscope ■ Store slides at RT

www.affymetrix.com/panomics

For research use only. Not for use in diagnostic procedures.