

QuantiGene ViewRNA miRNA ISH Cell Assay

24-Well Format for Suspension Cells

! **IMPORTANT:** If running QuantiGene ViewRNA miRNA ISH Cell Assay for the first time, please refer to the appropriate *QuantiGene ViewRNA miRNA ISH Cell Assay User Manual*, based on your sample type and processing format, to review assay guidelines and detailed procedures.

Part 1: Sample Preparation

Step	Task
Step 1. Coat Glass Coverslips with Poly-L-Lysine 1 hr	<ul style="list-style-type: none"> ■ Prepare 25 mL of 0.01% poly-L-lysine solution ■ Ensure availability of 100% ethanol, sterile 1X PBS, 24-well tissue culture plate, 12-mm glass coverslips, fine-tipped forceps ■ Sterilize coverslips in 24-well plate with 100% ethanol for 15 min at RT ■ Remove ethanol and air dry coverslips completely for 15 min at RT ■ Coat coverslips with diluted poly-L-lysine for 15 min at RT ■ Wash coverslips 3 times in 1X PBS ■ Air dry coated coverslips O/N at RT in laminar flow hood
Step 2. Culture Suspension Cells 20 min	<ul style="list-style-type: none"> ■ Culture cells under the recommended growth conditions until confluent

Part 2: Pretreatment and Target Probe(s) Hybridization

Step	Task
Step 3. Fix Samples 1 hr 15 min	<ul style="list-style-type: none"> ■ Prepare 1 L of 1X PBS, 4.5 mL of 4% formaldehyde in 1X PBS ■ Pellet cells, remove culture medium and wash once with 1X PBS ■ Fix cells in suspension with 4% formaldehyde for 1 hr at RT (in a fume hood), mix every 20 min ■ Pellet and wash cells 3 times with 1X PBS ■ Adjust cell concentration to 1000 cells/μL with fresh 1X PBS ■ Spot and evenly distribute 20,000 cells per poly-L-lysine coated glass coverslip ■ Bake at 50 °C for 30 min or until cells are completely dried
Step 4. Optional. Dehydrate/Rehydrate 10-20 min	<ul style="list-style-type: none"> ■ Ensure availability of 50%, 70%, and 100% ethanol if dehydrating cells for storage or for rehydrating cells for the assay ■ Dehydrate cells if storing sample: 100% EtOH and store at -20 °C ■ Rehydrate cells for assay if starting with dehydrated sample: 100% EtOH, 2 min -> 70% EtOH, 2 min -> 50% EtOH, 2 min -> 1X PBS, 10 min
Step 5. Prepare Buffers and Reagents during Sample Fixation 10 min	<ul style="list-style-type: none"> ■ Set a dry incubator to 40 ± 1 °C ■ Equilibrate EDC to RT ■ Pre-warm Probe Set Diluent QF to 40 °C ■ Thaw Probe Set(s) and place on ice ■ Place Protease QS on ice ■ Prepare 800 mL of Wash Buffer, 30 mL Storage Buffer ■ Ensure availability of Cross-Linking Buffer QM, Detergent Solution QC
Step 6. Cross-Link with EDC 1 hr 30 min	<ul style="list-style-type: none"> ■ Equilibrate cells with Cross-Linking Buffer QM 2 times, 10 min each at RT with rocking ■ Prepare 0.16M EDC just before use ■ Cross-link target with 0.16 M EDC for 1 hr at RT, rocking ■ Wash cells 3 times with 1X PBS

Step	Task
Step 7. <i>Optional.</i> Dehydrate/Rehydrate 10-20 min	<ul style="list-style-type: none"> ■ Dehydrate cells if storing sample: 50% EtOH 2 min -> 70% EtOH, 2 min -> 100% EtOH, 2 min -> Fresh 100% EtOH and store at -20 °C ■ Rehydrate cells for assay if starting with dehydrated sample: 100% EtOH, 2 min -> 70% EtOH, 2 min -> 50% EtOH, 2 min -> 1X PBS, 10 min
Step 8. Permeabilize and Digest with Protease 40 min	<ul style="list-style-type: none"> ■ Permeabilize cells with Detergent Solution QC for 10 min at RT ■ Wash cells 2 times with 1X PBS ■ Prepare Working Protease Solution ■ Treat cells with Working Protease Solution for 10 min at RT ■ Wash cells 3X with 1X PBS
Step 9. Hybridize Target Probe(s) 3 hrs	<ul style="list-style-type: none"> ■ Prepare Working Probe Set(s) Solutions ■ Hybridize cells with Working Probe Solutions for 3 hrs at 40 °C
Step 10. Wash Samples 10 min	<ul style="list-style-type: none"> ■ Wash cells 3 times in Wash Buffer, soaking 3 min each
Step 11. Stop Point 1 min	<ul style="list-style-type: none"> ■ Store samples in Storage Buffer, covered with lid and seal with parafilm at 4 °C. Do not exceed 24 hrs

Part 3. Signal Amplification and Detection

Step	Task
Step 12. Prepare Additional Buffers and Reagents 10 min	<ul style="list-style-type: none"> ■ Prewarm PreAmplifier Mix QM, Amplifier Mix QM, and Label Probe Diluent QF to 40 °C ■ Equilibrate the following reagents to RT: <ul style="list-style-type: none"> □ AP Enhancer, Fast Red Tablets, Naphthol Buffer □ Cell samples ■ Thaw, mix and place on ice 100X DAPI and Label Probe 4-Alexa 488 ■ Prepare 11 mL 1X DAPI and store on ice protected from light ■ Place Label Probe 1-AP on ice ■ Ensure availability of Wash Buffer and 1X PBS ■ Ensure availability of Prolong® Gold Antifade Reagent
Step 13. Wash Storage Buffer off Samples 5 min	<ul style="list-style-type: none"> ■ Wash cells 2X with Wash Buffer
Step 14. Hybridize with PreAmplifier 1 hr	<ul style="list-style-type: none"> ■ Prepare Working PreAmplifier Mix Solution ■ Hybridize cells with Working PreAmplifier Mix Solution for 1 hr at 40 °C
Step 15. Wash Samples 10 min	<ul style="list-style-type: none"> ■ Wash cells 3X with Wash Buffer, soaking 3 min each
Step 16. Hybridize with Amplifier 1 hr	<ul style="list-style-type: none"> ■ Prepare Working Amplifier Mix Solution ■ Hybridize cells with Working Amplifier Mix Solution for 1 hr at 40 °C
Step 17. Wash Samples 10 min	<ul style="list-style-type: none"> ■ Wash cells 3X with Wash Buffer, soaking 3 min each
Step 18. Hybridize with Label Probe(s) 1 hr	<ul style="list-style-type: none"> ■ Prepare Label Probe(s) Working Solution ■ Hybridize cells with Label Probe(s) Working Solution for 1 hr at 40 °C
Step 19. Wash Samples 20 min	<ul style="list-style-type: none"> ■ Wash cells 3 times in Wash Buffer, soaking 3 min each for the first 2 washes and 10 min for the final wash

Step	Task
Step 20. Fast Red Development 1 hr 20 min	<ul style="list-style-type: none"> ■ Incubate cells with AP-Enhancer for 5 min at RT; prepare Fast Red Substrate ■ Remove AP-Enhancer, add Fast Red Substrate and incubate at 40 °C for 45 min ■ Wash cells 2 times with 1X PBS ■ Fix cells with 4% formaldehyde for 10 min at RT (in a fume hood) ■ Wash cells 3 times with 1X PBS
Step 21. DAPI Counter Stain 5 min	<ul style="list-style-type: none"> ■ Stain cells with 1X DAPI for 5 min at RT ■ Wash cells once with 1X PBS and replace with fresh 1X PBS
Step 22. Mount Samples on Glass Slides 15 min	<ul style="list-style-type: none"> ■ Mount glass coverslips to microscope glass slide ■ Cure samples overnight at RT (protected from light) ■ View under microscope <ul style="list-style-type: none"> □ DAPI for nuclei: EX 387/11 nm, EM 447/60 nm (recommended Semrock P/N DAPI-11060B) □ ALEXA 488 for mRNA: EX 485/20 nm, EM 524/24 nm (recommended Semrock P/N FITC-2024B) □ Fast Red for miRNA: EX 530/40 nm, EM 590/40 nm (recommended Semrock P/N TRITC-B) □ ALEXA 750 for mRNA: EX 708/75 nm, EM 809/81 nm (recommended Semrock P/N Cy7-B)

www.affymetrix.com/panomics

For research use only. Not for use in diagnostic procedures.

P/N 19164 Rev. A 121004

© 2012 Affymetrix, Inc. All rights reserved. Affymetrix® and QuantiGene® are registered trademarks of Affymetrix, Inc. All other trademarks are the property of their respective owners.